

➤ Prefilled Syringe Manufacturing

Work with a Global Leader in Syringe Fill/Finish

Over **660M**
units manufactured
in the past 10 years



Global Experience

Experience manufacturing products sold worldwide including North America, South America, Europe and Asia.

Broad Expertise

Serving client needs with clinical through commercial manufacturing, customized support services, and experience with small molecules, biologics, and vaccines (adjuvant & conjugate).

Risk Mitigation Strategy

Like-in-kind high speed syringe lines located in different buildings with separate water and HVAC systems, creating a built-in risk mitigation strategy to help assure product availability.

Resources to Get You to Market

By providing high-speed, high-volume syringe filling, with automated inspection, labeling, and packaging equipment that mirrors filling capacity, BioPharma Solutions has the resources to get you to market.

Commitment to Our Clients

We Take Our Commitment Seriously

With over 60 pharmaceutical clients, BioPharma Solutions offers strong organizational capabilities to help ensure our collaboration provides the value you deserve.

Focusing on Results

The BioPharma Solutions team is committed to a common customer service philosophy - to provide service with a passionate focus on results.

Specialized Client Teams

BioPharma Solutions also staffs each project with a dedicated client team composed of cross-functional experts.



BioPharma Solutions' facility in Bloomington, Indiana, USA, is one of the largest contract manufacturers of prefilled syringes in North America.



↳ Prefilled Syringe Manufacturing

Work with a Global Leader in Syringe Fill/Finish

BioPharma Solutions' processes and systems, including high-speed contract manufacturing and automated inspection systems, are tailored to meet your parenteral manufacturing needs.

Syringe Capabilities

- Aseptic formulation and fill/finish
 - 0.5 mL - 20 mL syringe sizes
 - Diluents for Reconstitution
- Plastic or glass syringes that are sterile and ready-to-fill in a nested tub configuration, double-bagged
- Filling options: rotary piston pump or peristaltic filling mechanisms

Packaging Flexibility

- Option for manual or full automated Eisai inspection
- Option to choose from a variety of secondary packaging options and product enhancements including plunger rod, finger flange, product insert, carton, and blister
- Option to design custom label
- Serialization capabilities

Contact Us

Baxter Pharmaceutical Solutions LLC
927 South Curry Pike
Bloomington, Indiana, USA 47403
biopharmasolutions@baxter.com

Full-Service R&D Support

- Formulation development
- Silicone and/or tungsten sensitivity evaluation
- Container selection feasibility study & development stability studies
- Shear sensitivity and pumping mechanism studies
- Particulate matter measurement and characterization
- Extractables and leachables testing (outsourced)
- Functional testing (Glide force - break force)
- Analytical method development and validation

QC Stability

- Storage at ICH conditions in validated chambers
- Sample testing per stability protocol or project plan
- Assistance in stability project design
- Documentation through LIMS reports available

Regulatory

- Support the authoring of CMC documents for a client's submission
- Maintain and provide references to established Drug Master Files for:
 - Facility
 - Syringe container closure
- Additional support as needed



baxterbiopharmasolutions.com